

# Cambridge (CIE) A Level Chemistry



Your notes

## Infrared Spectroscopy

### Contents

\* Infrared Spectroscopy




## Interpreting IR Spectra

- **Infrared (IR) spectroscopy** is a technique used to identify compounds
  - It works by measuring how covalent bonds in molecules absorb infrared radiation
  - This absorption causes the bonds to vibrate in specific ways
- A **spectrophotometer** irradiates the sample with infrared waves
  - Then, it detects how much radiation is absorbed at each frequency
  - The resulting spectrum shows which types of bonds are present

## Bond vibrations and absorption

- All **covalent bonds** behave like tiny **springs** rather than rigid bars
- Like springs, they can **vibrate in different ways**, such as stretching and bending




**Different modes of vibration in molecules. Each mode has a characteristic frequency of vibration**

- These vibrations occur at **specific frequencies** in the **infrared (IR)** region of the electromagnetic spectrum
- When IR radiation matches a bond's natural frequency, the bond **absorbs energy** and vibrates more strongly
- Each **vibration mode** (e.g. symmetric stretch, bending) has its own **characteristic frequency**

- These absorbed frequencies are expressed as **wavenumbers** ( $\text{cm}^{-1}$ ), which are the **reciprocal of wavelength**

## Absorption patterns and functional groups

- Each organic compound has a **unique IR spectrum**
  - Particularly in the region below  $1500\text{ cm}^{-1}$ , known as the **fingerprint region**
- This region contains many small peaks caused by complex bond vibrations that are difficult to assign to specific groups
- However, comparing the fingerprint region to spectra from a database allows **exact identification** of a compound
- This is useful, for example, for identifying a particular isomer in a homologous series
  - They will all have the same functional groups but **different fingerprint regions**



*The superimposed spectra of propan-1-ol and propan-2-ol show the same bonds above  $1500\text{ cm}^{-1}$  but different fingerprint regions*

- Each bond type absorbs within a **characteristic wavenumber range**
- Absorptions vary in **width** (broad or sharp) and **intensity** (strong or weak)
- For example:
  - O-H bonds in alcohols and carboxylic acids give **broad** peaks due to hydrogen bonding
  - C=O bonds in carbonyl compounds give **sharp**, strong peaks.
- Comparing an unknown IR spectrum with data for known compounds helps identify **functional groups**



Your notes

## Absorption range of bonds table



Your notes

Bond	Functional groups containing the bond	Characteristic infrared absorption range (in wavenumbers) / $\text{cm}^{-1}$
C-O	Hydroxy, ester	1040 - 1300
C-C	Aromatic compound, alkene	1500 - 1680
C=O	Amide	1640 - 1690
	Carbonyl, carboxyl	1670 - 1740
	Ester	1710 - 1750
C $\equiv$ N	Nitrile	2200 - 2250
C-H	Alkane	2850 - 2950
N-H	Amine, amide	3300 - 3500
O-H	Carboxyl	2500 - 3000
	Hydroxyl	3200 - 3600

- Due to some absorption bands overlapping each other, other analytical techniques such as mass spectroscopy should be used alongside IR spectroscopy to identify an unknown compound





### Worked Example

Look at the two infrared spectra below and determine which one corresponds to propanone and which one to propan-2-ol.



Your notes



**Answer:**

- IR spectrum **A** is **propanone** and spectrum **B** is **propan-2-ol**.
- In IR spectrum **A** the presence of a strong, sharp absorption around  $1710\text{ cm}^{-1}$  corresponds to the characteristic  $\text{C}=\text{O}$ , carbonyl, group in a ketone.
- In spectrum **B** the presence of a strong, broad absorption around  $3200\text{--}3500\text{ cm}^{-1}$  suggests that there is an alcohol group present, which corresponds to the  $\text{--OH}$  group in propan-2-ol.